Ebola

Rapid lateral flow test

Test rapide

Schnelltest

Rapid lateral flow test for the qualitative detection of Ebola virus antigen VP 40 in serum and throat swab.
For Professional In Vitro Diagnostic Use Only

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Ⓞ GB Instruction for use

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Stada Diagnostik Ebola Rapid lateral flow test
Rapid lateral flow test for the qualitative detection of Ebola virus antigen VP 40 in serum and throat swab

For Professional In Vitro Diagnostic Use Only

REF 2.13.1.122

Intended use

Stada Diagnostik Ebola Rapid lateral flow test is a rapid visually read chromatographic immunoassay for the qualitative detection of Ebola virus antigen in human serum and throat swab.

It is intended to be used as an aid to detect Ebola hemorrhagic fever in case of infectious diseases with fever. It is intended to be used only together with clinical information/ methods and other in vitro diagnostic procedures.

Summary

Ebola hemorrhagic fever (Ebola HF) is a severe, often fatal disease in humans and nonhuman primates (monkeys, gorillas and chimpanzees) that has appeared sporadically. The virus is one of two members of a family of RNA viruses called the Filoviridae. There are four identified subtypes of Ebola virus. 90 % of the outbreaks are caused by Ebola Zaire.

Principle

Stada Diagnostik Ebola Rapid lateral flow test is a qualitative, rapid, membrane-based immunoassay for the detection of Ebola virus antigen VP 40 in human serum and throat swab.

After applying the sample to the sample well it can react with particles coated with anti-Ebola antibodies. The mixture migrates upward on the membrane by capillary action.

By reaching the test line region antigen bound to the particles can react with the capture reagent anti-Ebola antibody immobilized on the test line of the membrane and generate a colored line.

The presence of a colored line in the test region indicates a positive result (presence of Ebola antigen), while its absence indicates a negative result.

To serve as a procedural control, a colored line will always appear in the control line region indicating that reaction conditions and technology worked fine (proper volume of specimen has been added and membrane wicking has occurred).

Kit components/Content of the kit

Kit for 20 determinations

20 test cassettes in aluminium pouch including desiccant bag

The test contains anti-Ebola antibody-coated particles, anti-Ebola antibodies and anti-mouse IgG immobilized on nitrocellulose membrane.

20 transfer pipettes

20 sterile polyester tipped applicators/swabs

20 extraction tubes

1 work station for up to 8 extraction tubes

1 extraction buffer 10 ml, ready-to-use

The buffer solution contains detergent and preservative.

1 package insert

Equipment required but not provided

Timer, for measuring the time between application of the specimen and interpretation of the results after ten minutes.

Specimen collection device.

Storage and stability

Store Stada Diagnostik Ebola Rapid lateral flow test unopened at +4° C to +30° C (+39° F to +86° F) and use by the date given on the label.

Do not use any component past the labeled expiration date.

Do not freeze.

The test cassette is sensitive to light or humidity, therefore do not use test devices that have become wet or if the pouch is damaged or shows signs of damage!

After opening of the test cassette pouch the test has to be performed immediately, at least within 1 hour.

After opening, the buffer vial can be stored at +4° C to +30° C (+39° F to +86° F) for at least 3 months.

Specimen

Specimen should be collected, handled and shipped using appropriate infection control precautions for Ebola, other hemorrhagic fever viruses or other infectious diseases.

Shipping may be performed according to the policies of the shipping performer, customs regulations, and the requirements of the receiving laboratory. If specimens are to be shipped, they should be packed in compliance with federal regulations for transportation of etiologic agents.

Serum

Only use human serum collected and prepared according to the manufacturer's instructions for the specimen collection device.

Only clear, non-hemolyzed specimens may be used.

Serum may be stored at +2° C to +8° C (+35° F to +46° F) for up to 3 days. For long-term storage, specimens may be kept below -20° C/-4° F.

Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing.

Specimens should not be frozen and thawed repeatedly.

Throat swab

Use specimen at room temperature only (+15° C to +30° C/+59° F to +86° F). Only use fresh throat swab. Do not store. Do not freeze.

Warnings and precautions

Handle all specimens and used test components as though they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure. Appropriate biosafety practices must be used when handling specimens and test components. These precautions include, but are not limited to the following.

- Wear personal protective equipment: powder-free gloves, lab coat, eye protectors etc.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.
- Decontaminate and dispose all specimens, test components and other potentially contaminated materials in accordance with local regulations.

Testing procedure


Fig. 1: test cassette

Preparation

Allow test components and specimen to equilibrate to room temperature (+15° C to +30° C /+59° F to +86° F) prior to testing.

- Remove the test cassette from the sealed pouch and use it as soon as possible. The test should be performed within 1 hour after opening the pouch.

- Place the test cassette on a clean and level surface.
- Clearly indicate the patient/specimen identification on the test device to avoid confusions.
- Place the extraction tube in the work station and indicate the patient/specimen identification on the tube.

Throat swab

- Hold the extraction buffer bottle vertically and add **8 to 10 drops of the extraction buffer solution to the extraction tube**. Only use the buffer solution provided with this test kit.
- Insert the provided sterile swab into oral cavity.
- Very gently push the swab until resistance is met at level of the pharynx.
- Gently rotate the swab for a few times over the mucous membrane.
- Immediately insert the swab with collected throat swab specimen, compress the bottom of the tube and rotate the swab vigorously to mix the reagents for at least 1 minute. Press the swab against the inside of the tube, and expunge as much liquid as possible from the swab by pressing and rotating the fiber portion against the wall of the tube. Discard the swab.
- Let the mixture incubate for 1 to 15 minutes.
- Mix the contents of the tube again by gentle swirling.
- Now the mixture is ready to be tested. Fill the provided pipette with the mixture.
- Hold the pipette above the sample well of the test cassette and squeeze **6 drops into the well**. Wait until each drop is absorbed, before adding additional drops.
- Start the timer.
- Wait 10 minutes from addition of the specimen/buffer mixture and then read the results. Do not interpret the results after more than 15 minutes.

Serum

- Hold the extraction buffer bottle vertically and add **3 drops (approx. 100 µl) of the extraction buffer to the extraction tube**. Only use the buffer included in the test kit.
- Hold the provided transfer pipette vertically, transfer **3 drops of specimen (approx. 100 µl) into the extraction tube** and incubate for 1 to 15 minutes.
- Now the mixture is ready to be tested. Fill the provided transfer pipette with the mixture, hold the pipette vertically above the sample well of the test cassette and **drop the whole mixture into the sample well**. Wait until each drop is absorbed, before adding additional drops.
- Start the timer.
- Wait 10 minutes from addition of the specimen/buffer mixture and then read the results. Do not interpret the results after more than 15 minutes.

Quality control

A procedural control is included in the test. A colored line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

If the control line fails to appear (insufficient specimen volume or incorrect procedural techniques are the most likely reasons), the test result is INVALID (figure 4). Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Interpretation of results

POSITIVE: two distinct colored lines appear. One line should be in the control region (C) and another line should be in the test region (T) (figure 2).

Any visible colored line in the test region should be interpreted as positive, even if the line in the test region appears lighter or darker than the line in the control region.


Fig. 2: schematic illustration of a positive result

NEGATIVE: one colored line appears in the control region (C). No apparent line appears in the test region (T) (figure 3).


Fig. 3: schematic illustration of a negative result

INVALID: control line fails to appear (figure 4). If there is no colored line in the control region (C), even if a colored line appears in the test region (T), the result is invalid and the test should be repeated with a new test cassette.



Fig. 4: schematic illustration of invalid results

Limitations

Stada Diagnostik Ebola Rapid lateral flow test is designed for in vitro qualitative detection of Ebola antigen VP 40 in human serum and throat swab. Other body fluids and pooled specimens may not give accurate results, and should not be used.

Stada Diagnostik Ebola Rapid lateral flow test should only be used in samples from patients with suspected Ebola infections (diseases with fever) as an aid to detect Ebola hemorrhagic fever. It should not be used as the sole criteria for the diagnosis of Ebola infection.

Stada Diagnostik Ebola Rapid lateral flow test will indicate the presence of Ebola antigen concentrations $\geq 10^6$ PFU/ml in the specimen. Therefore it should not be used to exclude ebola infections.

Results should be interpreted by a physician only.

Results should be used only together with clinical information/symptoms/methods and other in vitro diagnostic procedures (e.g. PCR techniques).

The intensity of the colored line in the test region does not correlate to the titer of antigen in the specimen.

A negative result at any time does not exclude the possibility of Ebola infection.

A positive result should be confirmed using another method and the result should be evaluated together with the overall clinical evaluation.

False positive or false negative results may be produced by:

- insufficient specimen volume or incorrect procedural techniques
- using inappropriate specimen or materials (other buffer as contained in the test kit, other volumes of specimen and buffer as indicated)
- disregard of the incubation time of 1 minute of specimen and extraction buffer
- disregard of the waiting time of 10 minutes after applying the mixture to the test cassette

Low titers of Ebola might result in a faint line appearing in the test region after a prolonged time.

The “Hook effect” in case of very high levels of antigen in the specimen can lead to false negative results.

At present, evaluation of Stada Diagnostik Ebola Rapid lateral flow test is still in progress. Further data currently is being collected.

Performance characteristics

The performance of the Stada Diagnostik Ebola Rapid lateral flow test has been compared with a WHO-certified PCR test using African patient samples in Conakry, Guinea.

		STADA Diagnostik Ebola		total
		+	-	
PCR	+	12	1	13
	-	2	103	105
total		14	104	

Analytic sensitivity: \leq CT 23, Specificity: 98%

References

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do sub-clinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H Lübbbers B, Cararetto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immuno-sorbent assay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GPS-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Towner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwiya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Spécificité : 98 %

Références

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. J Virol Methods 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. Bull Soc Pathol Exot 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do sub-clinical filovirus infections occur worldwide? Med Microbiol Immunol 1992; 181:43-55

4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. Wkly Epidemiol Rec 2003; 78:285-9

5. Hartmann H Lübbers B, Cararetto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. J Immunol Methods 1993; 166:35-44

6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immuno-sorbent assay for Ebola virus antigens in tissues of infected primates. J Clin Microbiol 1992; 30:947-50

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Analytische Sensitivität: \leq CT 23

Spezifität: 98 %

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1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. J Virol Methods 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. Bull Soc Pathol Exot 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do sub-clinical filovirus infections occur worldwide? Med Microbiol Immunol 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. Wkly Epidemiol Rec 2003; 78:285-9
5. Hartmann H, Lübbers B, Carretto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. J Immunol Methods 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immuno-sorbent assay for Ebola virus antigens in tissues of infected primates. J Clin Microbiol 1992; 30:947-50

7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GPS-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42

8. Towner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwiya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. J Virol 2004; 78:4330-41

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Distributed by/Distribué par/Vertrieb durch: STADAPharm GmbH, 61118 Bad Vilbel, www.stada.de