




Bioäquivalenzstudie Candesartan STADA® 32 mg Tabletten vs. Atacand® PROTECT 32 mg Tabletten



	Candesartan STADA® 32 mg	Atacand® PROTECT 32 mg	90% Konfidenzintervall
maximale Plasmakonzentration (C _{max}) ng/ml:	251.95 ± 77.98	252.63 ± 105.13	94.59 – 113.17
Zeitpunkt der maximalen Plasmakonzentration (t _{max}) h:	3.79 ± 1.09	3.68 ± 0.99	
Fläche unter der Konzentrations- Zeit-Kurve (AUC _{0-t}) ng*h/ml:	2726.70 ± 565.43	2853.62 ± 742.54	92.03 – 101.85

**Candesartan STADA® 32 mg Tabletten ist bioäquivalent mit
Atacand® PROTECT 32 mg Tabletten.**